

### MICROPETTE MULTI CHANNEL



### **FEATURES:**

- 8 and 12 channel pipettes are appropriate for 96 well plates.
- Dispensing head rotates for effortless pipetting convenience.
- Individual piston and tip cone assemblies allowing easy repair and maintenance.
- Compound material made tip cone secures high sealing performance.
- Compatible with most universal tip brands.

TYPES	VOLUME
8 Channel Adjustable Volume Pipettes	0.5-10μL, 5-50μL and 50-300μl
12 Channel Adjustable Volume Pipettes	0.5-10μL, 5-50μL and 50-300μl

# **MICROPETTE (Mechanical)**





12 Channel Pipette

### **FEATURES:**

• Pipettes cover a volume range from 0.1µL to 10mL

**8 Channel Pipette** 

- Ergonomic design provides excellent operating experience
- Large display window allows for easy volume identification
- Easy calibration and maintenance
- Each MicroPette supplied with an individual calibration certificate according to ISO8655

# **MICROPETTE PLUS (Autoclavable)**



#### **FEATURES:**

- Fully autoclavable
- Ergonomic design provides excellent operating experience
- Easy-to-read volume display
- The pipettes cover the volume range from 0.1 µl to 5ml
- Easy calibration and maintenance
- Manufactured from innovative materials
- Each MicroPette Plus supplied with an individual calibration certificate according to ISO8655

TYPES	VOLUME
Autoclavable	0.1µl - 2.5µl, 0.5µl - 10µl, 2µl - 20µl, 20µl - 200µl, 5µl - 50µl, 20µl - 200µl, 50µl - 200µl, 10µl - 100µl, 100µl - 1000µl, 1000µl - 5000µl
Single Channel Fixed Volume	5րl, 10րl, 20րl, 50րl, 100րl, 500րl, 1000րl

#### **PIPETTOR STAND**

The linear and round stands have been designed to fit perfectly our pipette range, MicroPette and MicroPette Plus Autoclavable Pipette. The stands are convenient to hold up to 6 pipettes, both single and multi-channel.



## **ROUND**

Rotates smoothly



LINEAR

Static



