

The Orbital and Rocking Shakers provide a smooth continuous motion for uniform mixing. Display of time and speed can be viewed through the double LCDs. A wide range of platform is available for use with a variety of vessels.

SK-0180-S (Orbital)




FEATURES:


- Max. load capacity is 2 Kg. with platform
- Double LEDs for independent display of time and speed
- Speed range of 40-200rpm
- Wide range of platforms for use with a variety of vessels

SPECIFICATIONS	
Orbital Diameter	20mm
Max. Load Capacity (With Platform)	3Kg
Motor Type	DC Motor
Motor Input Power	16W
Motor Output Panel	10W
Power	220V, 50/60Hz
Time Setting Range	1 min-19h 59 min.
Speed Range	40-200rpm
Speed Display	LED
Timer Display	LED
Voltage	220V, 50/60Hz
Protection Class	IP21
Weight	3.1Kg
Permissible Ambient Temperature and Humidity	5-40°C 80% RH
Dimension [WxDxH]	270 x 330 x 110mm

SK-R1807-S (Rocking)







- Speed range of 10-80rpm
- Double LEDs for independent display of time and speed
- Wide range of platforms for use with a variety of vessels

SPECIFICATIONS	
Tilt Angle	7°
Max. Load Capacity (With Platform)	3Kg
Motor Type	DC Motor
Motor Input Power	16W
Motor Output Panel	10W
Power	220V, 50/60Hz
Time Settings Range	1min-19h 59min
Speed Range	10-80rpm
Speed Display	LED
Timer Display	LED
Voltage	220V, 50/60Hz
Protection Class	IP21
Weight	3.2 Kg.
Permissible Ambient Temperature and Humidity	5-40°C; 80% RH
Dimension [WxDxH]	270 x 330 x 130mm

