

Digital Centrifuge



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

This space-saving centrifuge provides reliable and consistent operation with its maintenance-free, brushless motor. The unit features two digital displays; one for speed and one for the timer.

Brand: DragonLab and Scilogex

Model: DM0412

Capacity: 12 Placer

Features:

- Speed range of 300-4500rpm
- User-friendly large LCD display provides all visual information
- Quiet and stable operation
- Brushless DC motor
- Dual protective casing

Reviews

There are yet no reviews for this product.