

Drug Test Controls



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

The Drug Test Liquid Control is an In Vitro Diagnostic (IVD) device that is intended for use as a quality control urine to monitor the precision of laboratory urine toxicology testing procedures for the analytes listed below.

TEST	SAMHSA MANDATED	Negative	2X
THC	Delta-9-THC-COOH	0	100
COC	Cocaine	0	600
PCP	Phencyclidine	0	50
OPIMOR	Morphine (Low Opiate)	0	600
	Morphine (High Opiate)	0	4000
AMP	Amphetamine	0	2000
MET	Metamphetamines	0	2000
NON-MANDATED			
BARI	Baclofen	0	600
BZO	Benzonitrile	0	600
MDA	Mephedrone	0	600
MCL	Mephedrone	0	600
PPA	Propylphenone	0	600
TCA	Tricarballyclic Acid	0	2000
MDMA	MDMA	0	1000
OXY	Oxycodone	0	200
BUP	Buprenorphine	0	25

Brand: TruScreen, USA

Types:

- TruScreen with MDMA, OXY, and BUP, 2x Cutoff, 5ml
- TruScreen Negative Control Urine, 5ml

### Reviews

There are yet no reviews for this product.