


## Fecal Occult Blood Rapid Card



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

### Description

A qualitative test for the detection of human blood hemoglobin in human stool specimens and early diagnosis for Colon Cancer.

**Brand:** Wondfo

### General Information:

- Screening test for Fecal Occult Blood
- Immunochromatography based one-step in vitro diagnostic test
- No instruments needed
- All materials provided, ready to use reagent
- No Cross reactivity - Bovine Hemoglobin, Swine hemoglobin, Horse hemoglobin, Bilirubin, Bovine serum albumin Vitamin C
- Analytical sensitivity : 50ng/ml of FOB
- Sensitivity : 99%(97/98), Specificity: 99.9%(202/202)
- Specimen : Stool
- Reading Time : 5~10 min
- Storage : 2~30°C

#### Reviews

There are yet no reviews for this product.