

Glucose Meter



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

A medical device for determining the approximate concentration of glucose in the blood.



Simple and Convenient

No coding required with GOD enzyme technology

Strip and Lancet ejectors for easy and safe disposal

Permits re-dosing of 2nd blood sample within 3 seconds

Backlight and strip port light for testing in dark settings

Quick and Reliable

Point of Care Testing : Glucose Meter

Accurate results in just 5 seconds with only 0.8 µL sample required

Operating temperature range of 5-45 °C allows for testing in almost any environment

Advanced Hct Auto-correction technology with a Hct testing range of 20 - 70%

Hypo, hyper, and ketone warnings remind user to take necessary action

Better Data Management

500 test memory with date and time

Easy PC Download using a USB cable

7, 14, 30, 60, and 90 day averaging allows user to quickly access trends in blood glucose levels

Meal markers included to monitor results before and after meals

510(k) cleared, CLIA Waived, CE Marked

Reviews

There are yet no reviews for this product.