

HCG



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

An immunochromatography based one-step in vitro diagnostic test of the qualitative detection of hCG in human urine.

Brand: Wondfo

General Information:

- Analytical sensitivity : 25mIU/ml of hCG
- Sensitivity : 99.9%, Specificity : 99.9%
- Specimen : Urine/Serum
- Reading Time : 3 min for Urine, 5 min for Serum
- Immunochromatography based one-step in vitro diagnostic test
- No cross reactivity: hLH(300mIU/mL), hFSH(300mIU/mL), and hTSH(1,000?IU/mL)
- Storage : 4~30°C

Reviews

There are yet no reviews for this product.