

Hemoglobin Meter



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Portable, low-cost hemoglobin meters that are user-friendly and reliable that measures hemoglobin levels in the blood.

Quick

Hemoglobin (Hb) and Hematocrit (Hct) results in < 15 seconds

Simply insert strip, apply specimen and read results

Accurate

Only 10 µL capillary or venous blood required

Precise results equivalent to automatic hematology analyzers

Wide Hb measurement range of 4.5-25.6 g/dL with Hct range of 13-75%

Stores up to 1,000 results

Portable

Point of Care Testing : Hemoglobin Meter

Operated by battery or optional AC adapter

Ready to use in any point of care setting

Dry strips eliminate inconvenience in shipping and storage of liquid reagents

Convenient

Quick data transfer via Mini USB port

Reviews

There are yet no reviews for this product.