

Histology Base Molds



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Metal Base Molds enable tissue to be positioned and embedded in wax when using routine or biopsy cassettes. These re-usable base molds are made from stainless steel and the internal rounded corners allow for easy removal of the formed wax block. They can be used on all standard size cassettes.

Brand: TruLab ®

Types:

- 7 x 7 x 7 mm
- 15 x 15 x 7 mm
- 24 x 24 x 7 mm
- 37 x 24 x 7 mm

Reviews

There are yet no reviews for this product.