


HIV 1/2



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

An in vitro diagnostic rapid chromatographic immunoassay for the qualitative detection of antibodies specific to HIV-1 including subtype O and HIV 2 in human serum, plasma or venous and capillary whole blood. It is intended as an aid to diagnosis of HIV in symptomatic, asymptomatic populations and persons at risk of HIV infection.

Brand: SD Biosensor, Korea

General Information:

- Immunochromatography based one-step in vitro diagnostic test
- No instruments needed
- Applied highly purified HIV recombinant proteins: p24, gp41 and gp36
- Detect HIV type 1 and type 2 with separate test lines at one time
- Sensitivity : 99.9%, Specificity : 99.7%
- Specimen : Serum, Plasma or Whole Blood
- Reading Time : 5~20 min
- Storage : 1~30°C
- WHO PQ Approved

Reviews

There are yet no reviews for this product.