

Immersion Oil



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

The BASO Immersion Oil is based on the formulation recommended by the World Health Organization for the following reasons: moderate viscosity, clear back ground, long lasting (will not dry out easily) and will not decolorize stained bacteria.

Brand: BASO (Taiwan)

Refractive Index: 1.5105~1.5125

Sizes:

4 x 20ml
30 ml
500 ml
1000 ml

Reviews

There are yet no reviews for this product.