

Incubator



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Lab incubators are necessary equipment for any laboratory conducting cell culture and tissue culture work. Incubators protect cells from changes in temperature, humidity, carbon dioxide and oxygen.

Brand: Digisystem

Features:

- Simple keypad input allows easy temperature setting
- LED digital display enables users to monitor the chamber temperature at any given moment
- Visual alarm indicator alerts users of abnormal conditions of the chamber temperature exceeds the setting points by 10°C
- Time setting can be continuous or timed mode (999 hours or 999 minutes)
- The temperature can be controlled and maintained to 70°C
- Shelves are adjustable
- Double doors allow sample viewing from the acrylic door

Reviews

There are yet no reviews for this product.

