

Innoculating Loop Sterilizer



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

DLab Infrared sterilizers are designed for sterilization of inoculating loops and needles.

Brand: Dlab

Model: ST800-S

- They are fast, convenient and can replace alcohol lamps and Bunsen burners
- Sterilization is completed within 5-7 seconds
- Very safe to operate with no visual fire
- LED display of actual temperature with temperature control accuracy of 10°C
- Can be used in biology safety cabinets and clean bench space
- High temperature accuracy with overheating protection
- Precise temperature control precision and adjustable from 300-830°C

Reviews

There are yet no reviews for this product.

