


Inoculation Loops



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Calibrated 1 and 10µL disposable loops, made of ABS material (soft flexible plastic), for a smooth, problem free planting and streaking of cultures. Our loops do not cut or gouge the agar surface during streaking.

Brand: TruLab ®

Size: Sterile, Individually Wrapped

- 1 µl + Needle, Neutral
- 10 µl + Needle, Blue
- 1 µl + 10 µl, Yellow
- Needle

Substitute: Inoculation Loop Holder with PVC Holder

Description: Loops are made of Nickel-Chrome metal wire. A strand of wire is molded around a specific pin size to form the loop shape and then the ends of wire strand are twisted together.

Complimentary Products:

- Nichrome Wire Loop
- Nichrome Wire Needle



Reviews

There are yet no reviews for this product.