

Insight U500 Urine Analyzer



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

The U500 Urine Analyzer is a semi-automatic and very easy to use with its touch screen feature. It can read up to 4 types of Urine Strips – 4, 5, 10 and 11 parameters.

Model: U500

Features:

Automated Features:

- 500 tests/hour allows efficient testing
- Automatic strip detection, alignment and waste disposal
- Automatically calibrates the analyzer whenever it is turned on.

Simple and Easy

- Simple menu accessible using large touch screen LCD
- Strip carrier/waste tray unit designed for one-step easy cleaning
- Routine, STAT and QC test modes available

Flexible Data Management

- Abnormal results flagged on print-out and highlighted on the screen automatically
- Memory capacity of 2,000 records and RS232C port for simple data transfer

Reviews

There are yet no reviews for this product.