

Insight Ultra U120 Urine Analyzer



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

The U120 Ultra Urine Analyzer is a high quality and competitively priced system developed for accurate and efficient urinalysis testing.

Model: U120

Features:

Easy to Operate

- Large color touchscreen LCD for simple menu navigation
- Up to 2,000 patient memory and 800 Operator ID storage using more than 20 types of strips
- Ability to select Time Logout between 1-99 with minutes or hours option

Accurate and Efficient

- Advanced CMOS Image Sensor ensures accurate readings
- Can read strips with up to 11 parameters, including Macroalbumin/Creatinine
- Ability to edit abnormal results

Simple Data Transfer

- Bi-directional transmission of LIS data using Bluetooth, LAN or WLAN
- Ability to update software with SD card or USB flash drive

Unique Lockout Functions

- Strip, User and QC Lockout

Reviews

There are yet no reviews for this product.