

Liquid Urine Controls



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Simply validate visual and urinalysis with Insight Liquid Urine Controls.

Features:

Reliable

- Use with Mission and Mission Expert Urinalysis Reagent Strips and Urine zers for optimum quality control
- Validate urinalysis results and prevent procedure errors
- Control Level 1 provides negative results for LEU, NIT, URO, PRO, pH, BLO, SG, KET, BIL, GLU and ASC
- Control Level 2 provides positive results for LEU, NIT, URO, PRO, pH, BLO, SG, KET, BIL and GLU with negative results for ASC

Quick and Convenient testing

- Ensures accurate results for all parameters
- Obtain quick results in any setting
- Competitively priced

Packaging and Shelf Life

- Ready-to-use without dissolving in distilled water
- 18 months shelf life for unopened controls at 2-8 C
- Dropper Tip Bottles - available in separate positive and negative levels
- Dropper tip bottles provide efficient use of the control solution
- Easily drop the control solution onto each reagent pad using the dropper tip bottle
- Control can be used up to 40 times within 20 days at room temperature

Reviews

There are yet no reviews for this product.