

MET/THC



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

MET/THC test is a rapid, one step, immunochromatographic assay designed for qualitative detection of drug metabolite in human urine at a cut-off concentration.



Brand: Wondfo

General Information:

- Immunochromatography based one-step in vitro diagnostic test
- Conformed to USA, European, and International drug of abuse test standards
- No instruments needed
- Suitable for Point-of-Care Testing (POCT)
- Cut-off: 50ng/ml

Drug Testing Kits : MET/THC

- Sensitivity : 100%
- Specificity : 100%
- Specimen : Urine
- Reading Time : 5 min
- Storage : 4-30°C

Reviews

There are yet no reviews for this product.