

## Mounting Medium



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

### Description

The Eukitt mounting medium (contains 45% acrylic resin and 55% xylenes), used in microscopy. It is an adhesive and specimen preservative that can be used manually and in automated coverslipping equipment. Special purpose adhesive and sealant for use with glass, non plastic and other non-reactive materials. It spreads quickly and evenly without forming any air bubbles. Stays optically clear for over 10 years.

**Brand:** Eukitt (Germany)

**Refractive Index:** 1.5105

**Sizes:**

500 ml

1000 ml

### Reviews

There are yet no reviews for this product.

