

Paraffin Wax Trimmer



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Wax Trimmer is an effective way of removing excess paraffin wax from the tissue cassettes after embedding.

Features:

- High quality heating element, delivers prompt heating and energy saving.
- US made temperature sensor enables the unit work with high precision and reliability
- Auto-memory of last settings and recovery functions
- LED display of heating status
- Intelligent heating control of wax valve ensures smooth wax flow

Reviews

There are yet no reviews for this product.