

Pipette Shaker



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

The Pipette Shaker can easily achieve the uniform mixing of the chemicals with blood while they are in the diluting pipettes.

Brand: Digisystem

Model: PS-600T

Features:

- Has a pipette holder made of stainless steel spring clips.
- Provides continuous or timed operation (0-15 minutes)
- Can hold up to 6 pipettes

Reviews

There are yet no reviews for this product.