

Sharps Container Systems



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Sharps container systems provide solutions for the safe and effective containment and disposal of contaminated needles, syringes, and other sharps. Containers are stackable, puncture-resistant, and autoclave- and incinerator-safe. They feature transparent, locking lids and needle removal ports. The tortuous path systems restrict hand entry.

Brand: TruLab ®

Types: Containers

1L -



5L -



Needle Destroyers

- With cutter, Round
- With cutter, Square



Reviews

There are yet no reviews for this product.