

Slide Dispenser



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Easy-to-use top-loading dispenser accepts standard 25x75mm (1x3") slides. Turning either right or left side wheel dispenses slides individually while keeping the rest of the slides dust free.

Brand: TruLab ®

Capacity: 91 Slides

Reviews

There are yet no reviews for this product.