Slide Staining Plate



Rating: Not Rated Yet Price:
Variant price modifier:
Base price with tax:
Price with discount:
Salesprice with discount:
Sales price:
Sales price without tax:
Discount:

Tax amount:

Ask a question about this product

Description

Ideal for Hematology, Cytology, and Microbiology labs, this staining plate has rails covered with a special "tacky" polymer strip that hold 20 slides in place, even if the tray is held at an angle. The raised rails facilitate easy retrieval of slides while the 4 rubber feet ensures greater stability. Comes with a cover made of PET to enable immediate visual observation.



Reviews

There are yet no reviews for this product.