


### Slide Staining System



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

#### Description

The Slide Staining System comes with 2 components: Staining Jar and Rack. The Staining Jar is a high-quality glass jar commonly used for all routine and special staining procedures. Each jar comes with a separate fitted cover which prevents evaporation and contamination and helps contain chemical odors. The side walls are slightly tapered to facilitate use with staining racks.

- 30 Slide Unit
- 20 Slide Unit

#### Reviews

There are yet no reviews for this product.