Specimen Seal up Racks



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

Ask a question about this product

Description

Seal Up Rack: A safe solution for blood specimen post processing.

Features:

Isolated Seal

» Avoids cross infection, ensuring safety of specimen, personnel and environment.

Multifunctional

» It can be separated into two parts. The bottom part can be used as universal tube rack in laboratory. No need to transfer repeatedly.

Number Marks

» Identify specimen easily with numbers (1~100) marks on the rack.

One Step Operation

» Seal up specimen in one step, improving work efficiency.

Stackable

» Positioning groove and ridge design allows stable mounting of racks on top of each other.

Reviews

There are yet no reviews for this product.