

Standard Q Antigen Test Kit Nasal Swab



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

Brand: Standard Q Antigen Test Kits Nasal Swab (SD Biosensor)

Features:

- Fast Results within 5 minutes
- Easy to use
- Specimen: Nasal Swab
- Sensitivity: 97.12%, Specificty: 100%
- All necessary reagents provided and no euipment needed
- FDA Approved

Performance Characteristics

> STANDARD Q COVID-19 Ag Test (Nasal) result by Central Research Laboratory in India

Type of specimen		RT-PCR		
		Positive	Negative	Total
STANDARD Q COVID-19 Ag Test (Nasal)	Positive	101	0	101
	Negative	3	399	402
	Total	104	399	503
Sensitivity: 97.12%, Specificity: 100%				

Reviews

There are yet no reviews for this product.