


Syphilis



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

A qualitative screening test kit for the detection of antibodies to *Treponema pallidum* in human blood.

Brand: SD Biosensor, Korea

General Information:

- *Treponema pallidum* Antibody Test
- Immunochromatography based one-step in vitro diagnostic test
- No instruments needed
- Sensitivity : 100% (169/170)
- Specificity : 100% (199/200)

Rapid Test : Syphilis

- Specimen : Serum, Plasma or Whole blood
- Reading Time : 5~20 min
- Storage : 1~30°C

Reviews

There are yet no reviews for this product.