

## TB



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

### Description

Intended to detect Mycobacterium tuberculosis (TB) MPT64 antigen in solid culture or liquid culture samples using immunochromatography.

**Brand:** SD Biosensor, Korea

### General Information:

- Immunochromatography based one-step in vitro diagnostic test
- No instruments needed
- Sensitivity : 100% (66/80)
- Specificity : 100% (306/330)
- Specimen : Solid Culture, Liquid Culture
- Reading Time : 10~15 min
- Storage : 2~30°C

### Reviews

There are yet no reviews for this product.

