

Test Tube Caps



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

These test tube caps are designed for use with 12 and 16mm diameter test tubes. The flanged threads allow the plug caps to descend into tubes until a friction fit is secure and plug is stationary.

Colors:

- Purple
- White

Sizes:

- For 5ml test tubes
- For 10ml test tubes

Reviews

There are yet no reviews for this product.