


## Troponin I



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

### Description

Intended for the qualitative determination of Cardiac Troponin I

: Wondfo

Brand

:

General Information

- Immunochromatography based one-step in vitro diagnostic test
- Point of care test at the Emergency Room
- No instruments needed
- No cross reactivity
- Specimen : Serum, Plasma or Whole Blood
- Reading Time : 15~20 min
- Storage : 2~30°C

	Troponin I	CK-MB	Myoglobin
Analytical Sensitivity	1ng/ml	5ng/ml	75ng/ml
Sensitivity	99.9%	99.9%	99.9%
Specificity	98.1%	94.6%	93.6%

#### Reviews

There are yet no reviews for this product.