Urine and Stool Collection



Rating:	Not Rated	Yet
Drice:		

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

Ask a question about this product

Description

TruLab ®oratories offers an extensive range of urine and stool containers. The containers are made entirely of polypropylene and polysterene, are chemical-resistant even to hydrochloric acid, and are metal-, latex-, zinc-, and fluorescence-free. Containers can be subjected to freezing, thawing, and EtO gas sterilization without causing changes in materials or physical appearance.

Brand: TruLab ®

Accessory: Specimen Container Rack

Types:

- Sterile
- Non Sterile
- Volume Capacity: 30, 60, 90, and 150 ml



This 24 placer, stainless steel rack keeps specimen containers organized and in place, minimizing risks of spilling and contamination.

Reviews

There are yet no reviews for this product.