

## Vortex



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

### Description

The Vortex is a device used in laboratories to mix liquid components in tubes and resuspend cells. It is widely used in biological and chemical analysis applications.

**Brand:** Dlab

**Model:** MX-S

### Features:

- Touch operation or continuous mode
- Variable speed control from 0-2500 rpm
- Used for various mixing applications with optional adapters
- Specially designed vacuum suction feet for body stability
- Robust aluminum-cast construction

#### Reviews

There are yet no reviews for this product.