

Wash Bottle



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Description

These bottles save time, help avoid mistakes, and increase safety. Bottles are durable, chemical-resistant, squeezable, and easy to clean. The leak proof wash cap assembly will deliver a single drop or a needle thin stream, with finger-tip squeeze control of flow. For a coarser flow, just cut the nozzle tip back.

Brand: TruLab ®

Sizes:

- 250 ml
- 500 ml

Reviews

There are yet no reviews for this product.

