



DENGUE RAPID TEST



DENGUE IgG/IgM RAPID TEST

Dengue IgG/IgM Rapid Test is a rapid chromatographic immunoassay for the qualitative and differential detection of IgG and IgM antibodies to Dengue virus in human serum, plasma or whole blood.

- Detection: IgG and IgM antibodies to Dengue virus
- Specimen: Serum/Plasma/Whole Blood
- Sensitivity: 96.7%
- Specificity: 98.6%
- Test Time: 15 minutes
- Shelf life and storage temperature: 24 months at 2-30°C

DENGUE NS1 ANTIGEN RAPID TEST

Dengue NS1 Antigen Rapid Test is a rapid chromatographic immunoassay for the qualitative detection of Dengue virus NS1 antigen in human serum, plasma or whole blood.

- Detection: Dengue NS1 Ag
- Specimen: Serum/Plasma/Whole Blood
- Sensitivity: 99.1%
- Specificity: 99.5%
- Test Time: 15 minutes
- Shelf life and storage temperature: 24 months at 2-30°C


DENGUE NS1 Ag & IgG/IgM COMBO RAPID TEST

Dengue NS1 Ag & IgG/IgM Combo Rapid Test is a rapid chromatographic immunoassay for the qualitative detection of Dengue virus NS1 antigen and IgG/IgM antibodies to Dengue virus in human serum, plasma, or whole blood.


- Detection: Dengue NS1 Ag and IgG/IgM antibodies to Dengue virus
- Specimen: Serum/Plasma/Whole Blood
- Test Time: 15 minutes
- Shelf life and storage temperature: 24 months at 2-30°C

TEST PROCEDURES

Dengue NS1 Ag



Dengue IgG/IgM



Dengue NS1 Ag & IgG/IgM Combo

