



SALIVA ALCOHOL TEST STRIPS

Mission® Saliva Alcohol Test Strips, a rapid saliva test for the determination of relative Blood Alcohol Concentration (BAC). Results in 2 minutes.



BROAD AND COMPLETE DETECTION

- Single strip detects multiple concentration of 0.020-0.30% BAC
- Unique strip with all-in-one detection offers cost-effective testing of a single sample
- 1 year shelf life

EASY TO USE

- Non-invasive specimen collection allows for on site testing
- No additional equipment or training necessary
- Easy to interpret results in 2 minutes.

ACCURATE AND RELIABLE

- Quick and easy detection
- Suitable for use in homes, schools, workplaces, rehabilitation centers, hospitals, military, law enforcement units and other government agencies.

CONVENIENT

- Cost-effective method compared to other alcohol detection test.
- Compact and individually wrapped packaging for use anytime and anywhere.

FEATURES	SPECIFICATIONS
Product Name	Saliva Alcohol Test Strips
Blow Bag	N/A
Technology	Enzymatic reaction
Specimen	Saliva
Blood Alcohol Concentration (BAC)	0.02%-0.30%
Measuring Units	% BAC
Reading Time	2 Minutes
Reading Stability	3 minutes
Storage temperature	2-27°C (36-86°F)
Shelf Life	1 year
Pouch Dimension (L x W)	13.5cm x 4.7 cm (5.3" x 1.9")
Weight	225.4 g (80 oz.)
Kit Box Dimension	18 cm x 10.5 cm x 5 cm (7.4" x 4.1" x 2.0")
Package Size (Test)	25
Carton Dimensions (L x W x H)	55 cm x 39.5 cm x 46 cm (22" x 16" x 18")
Weight	9.6 kg. (21.2 lbs.)
Number of Kits per Carton	80