

URINE STRIPS

The Insight® Urinalysis Strips are a great tool to get insight into the patient's overall health and their metabolic, liver and kidney functions. It is a very cost effective screening tool. The Mission® Urinalysis Strips are very reliable and have comparable performance to the market leaders.



AVAILABLE IN:

- 4 Parameters
- 4 parameters, Xpert
- 10 Parameters
- 10 Parameters, Xpert
- 11 Parameters

EASE OF USE

- Compatible for visual and analyzer reading
- Up to 11 combinations available for visual reading
- Up to 7 combinations available for analyzer reading

RELIABLE RESULTS

- High ascorbic acid resistance
- Prevents false-negative results for Glucose, Blood, Nitrite and Bilirubin
- Equivalent to the market leader for enhanced accuracy
- Analytical sensitivity better than or equivalent to market leaders
- High quality color chart ensures accurate visual reading
- Results in 2 minutes

CONVENIENT

- 100 strips per kit
- 2 years shelf life for closed canisters and 3 month shelf life for opened canisters

RAPID ANALYSIS OF URINE SPECIMENS:

- LEU
- PRO
- SG
- KETONE
- NIT
- pH
- GLU
- URO
- BLO
- BIL

LIQUID URINE CONTROL LEVEL 1 AND 2

Simply validate visual and urinalysis with Insight Liquid Urine Controls.



RELIABLE

- Use with Mission® and Mission® Expert Urinalysis Reagent Strips and Urine zers for optimum quality control
- Validate urinalysis results and prevent procedure errors results for ASC

QUICK AND CONVENIENT TESTING

- Dropper Tip
- Dropper tip bottles provide efficient use of the control solution
- Easily drop the control solution onto each reagent pad using the dropper tip bottle.

PACKAGING AND SHELF LIFE

- Ready-to-use without dissolving in distilled water
- 18 months shelf life for unopened controls at 2-8 °C
- Dropper Tip Bottles - available in separate positive & negative levels
- Dropper tip bottles provide efficient use of the control solution
- Control can be used up to 40 times within 20 days at room temperature

FEATURES	SPECIFICATIONS
Product Name	Liquid Urine Control
Test Parameters	LEU, NIT, URO, PRO, pH, BLO, SG, KET, BIL, GLU, ASC, ALB, CRE, CA (13)
Solution Detection Levels, Level 1	Negative: LEU, NIT, URO, PRO, pH, BLO, SG, KET, BIL, GLU, ASC, ALB, CRE, CA
Solution Detection, Levels, Level 2	Positive: LEU, NIT, URO, PRO, pH, BLO, SG, KET, BIL, GLU, ALB, CREA, CA; and Negative ASC
Compatible Urine Strips	Mission® Urinalysis Reagent Strips, Mission® Expert Urinalysis reagent strips
Reading Time / Stability	Refer to insert
Storage Temperature	2-8°C
Memory	24 Months
Opened Control Stability	30 days at 15-30°C or until the expiration date at 2-8°C
Maximum Test per Unit	20 or 40 tests/bottle