



DIGITAL CENTRIFUGE DM0424

FEATURES:

- Max. speed 4000rpm.
- Rotor speed can be set in RPM/RCF, with 3 programs storage capacity for routine applications.
- Max. capacity 50ml x 4, Adapt to different sizes of blood collection tubes and centrifuge tubes suitable for swing out rotor and angle rotor.
- Automatic internal diagnosis and self-detection.
- Maintenance-free brushless DC motor.



SPECIFICATIONS	DM-0424
Speed Range	500-4000 rpm Increment: 100 rpm
Max. RCF	2500 x g
Speed Accuracy	±100 rpm
Rotor	Swing 6 x 15ml Swing 4 x 50 ml Angle 24 x 15 ml Angle 12 x 10 ml / 8 x 15 ml
Run Time	1-99 min, HOLD
Motor	Brushless DC Motor
Display	LCD
Program	Blood: 3200 rpm, 10 mins. Urine: 1800 rpm, 5 mins. Faecal: 1300 rpm, 10 mins.
Safety Device	Automatic Internal Diagnosis
Acceleration / Deceleration	3↑4↓steps
Power	Single Phase 240V, 50Hz/60Hz
Dimension (W x D x H)	364 x 440 x 268 mm
Weight	14.5 kg
Noise	≤62dB(A)



TruLaboratories
Corporation

+632 82-TRULAB (878522)
+632 8287.3683
+63917.878.2288 • +63925.878.2288
www.trulaboratories.com