



DIGITAL CLINICAL CENTRIFUGE DM0408

FEATURES:

- Speed range of 300-4000rpm.
- 3 programs storage capacity for routine applications.
- Automatic internal diagnosis.
- Maintenance free brushless DC motor.



SPECIFICATIONS	DM-0408
Speed Range	300-4000rpm : 100 rpm
Max. RCF	1900 x g
Speed Accuracy	±100rpm
Rotor	Angle 8x15ml (Default) Angle 10ml x 12
Run Time	1-99 min, HOLD
Motor	Brushless DC Motor
Display	LED
Program	Blood: 3200 rpm, 10 min. Urine: 1800 rpm, 5 min. Faecal: 1300 rpm, 10 min.
Safety Device	Automatic Internal Diagnosis
Acceleration / Deceleration	25s↑35s↓
Power	Single-Phase 220V-240V, 50Hz/60Hz
Dimension (W x D x H)	286 mm x 367 mm x 227 mm
Weight	9.2 kg.
Noise	≤60dB(A)



TruLaboratories
Corporation

+632 82-TRULAB (878522)
+632 8287.3683
+63917.878.2288 • +63925.878.2288
www.trulaboratories.com