

BaSO

FECAL OCCULT BLOOD TEST SEMI QUANTITATIVE METHOD

A Fecal occult Blood Test is a non-invasive test (nothing enters the body). This test detects hidden (occult) blood in the stool. Such blood may come from anywhere along the digestive tract. Hidden blood in stool is often the first, and in many cases the only, warning sign that a person has a colorectal disease, including colon cancer.



Contents:

BaSO Fecal OB-II

1. Fecal OB-II Test Card
2. Developer A
3. Developer B
4. Positive Control
5. Applicator Stick

Principle :
Fecal OB-T Reagent kit

(Modified Pyramidon method)

- "Two function in One". Combine Fecal specimens collection tube and Test tube function. Cost savings.
- "Simple steps to collect specimens". Patient self-collected specimens will not have too much stool specimens.
- "Easy test and read test result in two minutes". Reaction in closed tube, clean and sanitary no bad odor.
- "Semi-quantitative and no carcinogen". Using the test results color change to achieve semi-quantitative purposes. Pyramidon method is safe for operator.

The Fecal OB-T reagent kit is a "Two in One" test kit, application on Fecal Occult Blood test.

Procedure:

