

DLAB

VORTEX MIXER

MX-S

FEATURES:

- Touch operation or continuous mode
- Variable speed control from 0 to 2500rpm (MX-S)
- Used for various mixing applications with optional adapters (MX-S)
- Specially designed vacuum suction feet for body stability
- Robust aluminum-cast construction



SPECIFICATIONS	MX-S (ADJUSTABLE SPEED)
Voltage	220V, 50/60Hz
Power	60W
Mixing Motion	Orbital
Orbital Diameter	4mm
Motor Type	Shaded Pole Motor
Motor Rating Input	58W
Motor Rating Output	10W
Speed Range	0-2500rpm
Speed Display	Scale
Run Type	Touch operation / Continuous
Dimension [WxHxD]	127 x 130 x 160mm
Weight	3.5kg
Permissible Ambient Temperature	5-40°C
Permissible Relative Humidity	80%RH
Potection Class	IP21

Accessories (Optional)



Standard Top (Default)
for <Ø30mm tubes and small vessels



Tube Holding Rod
used with tube adapters



Universal Top Plate
Ø100mm



Platform Pad
<Ø99mm tubes and small vessels



Vacuum Chuck
Made of Rubber

TUBE ADAPTERS

REMARKS



48 holes test tubes
Ø6mm



15 holes test tubes
Ø10mm



12 holes test tubes,
Ø12mm



for 8 holes test tubes
Ø16mm



for 8 holes test tubes
Ø20mm