
























## BLOOD COLLECTION TUBE

Blood Collection Tubes comes with various sizes, and color coded tops indicating tube contents. Most blood collection tubes contain an additive that either accelerates clotting of the blood (clot activator) or prevents blood from clotting (anticoagulant).






Mini Blood Collection Tubes are ideal for collecting and transporting capillary blood. It is highly recommended in patients with fragile veins, young children, geriatrics and in situations where only a small blood sample is required.

PICTURE	COLOR	ADDITIVE	ROLE OF ADDITIVE
	Purple / Lavender	EDTA (Ethylenediaminetetraacetic Acid)	Binds and remove calcium to prevent blood from clotting
	Light Blue	Sodium Citrate /9NC	Binds and remove calcium to prevent blood from clotting
	Green	Lithium Heparin Sodium Heparin	Inhibits thrombin formation to prevent clotting
	Red	Plain Tube or contains silica particles which act as clot activators	Clot activator promotes blood clotting with glass or silica particles
	Yellow	Anticoagulant SPS (Sodium Polyanetholsulfonate) and ACD (acid citrate dextrose)	Prevents the blood from clotting and stabilizes bacterial growth
	Black	Sodium Citrate	Forms calcium salts that remove calcium and inhibits clotting






COLOR	DESCRIPTION	PACKAGING
	EDTA, K2, 2ml, PS	100 pcs Tubes per Box
	EDTA, K3, Glass, 2ml	100 pcs Tubes per Box
	EDTA, K2, PS, 3ml	100 pcs Tubes per Box
	EDTA, K3, Glass, 3 ml	100 pcs Tubes per Box
	EDTA, K2, PS, 5ml, Short and Tall	100 pcs Tubes per Box
	EDTA, K3, Glass, 5ml, Short and Tall	100 pcs Tubes per Box
	Plain, Glass, 3ml	100 pcs Tubes per Box
	Plain, PS, 5ml, Short and Tall	100 pcs Tubes per Box
	Plain, Glass, 5ml, Short and Tall	100 pcs Tubes per Box
	Glass, 3.2% Sodium Citrate, 1.8ml	100 pcs Tubes per Box
	PS, 3.2% Sodium Citrate, 1.8ml	100 pcs Tubes per Box
	Lithium Heparin, Glass, 4ml	100 pcs Tubes per Box
	Sodium Heparin, Glass, 4ml	100 pcs Tubes per Box
	ESR, Glass, 3 ml	100 pcs Tubes per Box
	Glass, 4ml	100 pcs Tubes per Box
	PS, 4ml	100 pcs Tubes per Box
	PS, 5ml, 13 x 75mm	100 pcs Tubes per Box



COLOR	DESCRIPTION	PACKAGING
	PP Cap, Lithium Heparin, 0.5ml	50 pcs Tubes per Box
	PP and Rubber Cap, EDTA, K2 and K3, 0.5ml	50 pcs Tubes per Box
	PP Cap, 0.5ml	50 pcs Tubes per Box
	PP Cap, Gel and Clot Activator, 0.5ml	50 pcs Tubes per Box
	PP Cap, Gel and Clot Activator, Amber Tube, 0.5ml	50 pcs Tubes per Box



 +632 82-TRULAB (878522)  
 +632 8287.3683  
 +63917.878.2288 • +63925.878.2288  
 www.trulaboratories.com

## BLOOD COLLECTION TUBE



COLOR	DESCRIPTION	PACKAGING
	3.2% Sodium Citrate, Glass, 1.8ml	100 pcs Tubes per Box
	3.2% Sodium Citrate, PS, 1.8ml	100 pcs Tubes per Box
	Glass, 10 ml	50 pcs Tubes per Box
	Lithium Heparin, Glass, 4ml	100 pcs Tubes per Box
	Sodium Heparin, Glass, 4ml	100 pcs Tubes per Box
	Glass, 2.4ml	100 pcs Tubes per Box
	Glass, 4ml	100 pcs Tubes per Box
	PS, 4ml	100 pcs Tubes per Box
	Glass, 5ml, Tall Tube	100 pcs Tubes per Box
	PS, 5ml, Tall Tube	100 pcs Tubes per Box
	EDTA, K2, PS, 2ml	100 pcs Tubes per Box
	EDTA, K3, Glass, 2ml,	100 pcs Tubes per Box
	EDTA, K2, PS, 3ml	100 pcs Tubes per Box
	EDTA, Glass, K3, 3ml	100 pcs Tubes per Box
	EDTA, PS, K2, 5ml, Short	100 pcs Tubes per Box
	EDTA, Glass, K3, 5ml, Short	100 pcs Tubes per Box
	Plain, Glass, 3ml	100 pcs Tubes per Box
	Glass, with Procoagulation, 3ml	100 pcs Tubes per Box
	Plain, Glass, 4ml	100 pcs Tubes per Box
	Plain, PS, 4ml	100 pcs Tubes per Box
	Plain, Glass, 5ml, Short	100 pcs Tubes per Box
	Plain, PS, 5ml, Short	100 pcs Tubes per Box
	Glass, with Procoagulation, 5ml, Short Tube	100 pcs Tubes per Box
	Plain, Glass, 10ml	100 pcs Tubes per Box



COLOR	DESCRIPTION	PACKAGING
	EDTA, K2, Glass, 2 ml	100 pcs Tubes per Box
	EDTA, K3, Glass, 3 ml	100 pcs Tubes per Box
	EDTA, K3, Glass, 5ml, Tall Tube	100 pcs Tubes per Box
	PP Cap, 0.5 ml	100 pcs Tubes per Box
	Glass, 5 ml, Short Tube	100 pcs Tubes per Box
	3.2% Sodium Citrate, Glass, 2 ml	100 pcs Tubes per Box
	Plain, Glass, 4 ml	100 pcs Tubes per Box