



### EMBEDDING RING

These rings hold and identify tissue specimen blocks during sectioning, storing, positioning and embedding operations. In addition to this the provided embedding ring ensures accurate positioning of tissue samples.



### HISTOLOGY BASE MOLDS

For use in the preparation of pathology specimens. Available in stainless and plastic.

#### Sizes Available:

- 7 mm x 7 mm x 7 mm
- 15 mm x 15 mm x 7 mm
- 24 mm x 24 mm x 7 mm
- 37 mm x 24 mm x 7 mm



### FILTER PAPER

Filter paper are highly useful in solvents filtering cellulose analysis, and are a must have to any lab or industrial process that involved separation and purification.

Format: Ashless

Weight: 80 grams per sq. meter

Size 600 mm x 600 mm




### TISSUE AND BIOPSY CASSETTE

Disposable plastic cassettes hold tissue specimen during the embedding process, as well as in a storage life. Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. Available Colors in Yellow, Blue, Green, Pink and White.



### TALLY COUNTER

This hand tally counter records everything from 0000 to 9999, it is chrome plated case, lightweight with finger ring for the best taken in hand.



### MICROPLATES

Trulab®oratories offer a variety of 96 well microplates for diagnostic & immunological research applications. Available in U Bottom, V Bottom and Flat Bottom, all items are Sterile.