

BLUE PIPETTE TIPS

PIPETTE TIP

500 pcs. per pack

TYPES:



Eppendorf Gilsons with Crown with Crown with Crown



Gilsons

YELLOW PIPETTE TIPS



Eppendorf

TYPES:





with Crown

without without Crown with Graduation Crown

1000 pcs. per pack

WHITE PIPETTE TIPS



TYPES:



With Crown $0.5 - 10 \mu L$



With Crown 5000 µL

 $1000 \text{ pcs. per pack } (0.5 - 10 \mu\text{L})$ 100 pcs. per pack (5000 μL)



100 pcs. per pack

BRAND: DUB TYPES:



CONDUCTIVE TIPS

Conductive tips are produced under clean room conditions, they are tested to be from Dnase, Rnase and pyrogens. With conductive tips, you can use the smallest amount of liquid and get the highest precision without any contamination.

- 96 tips in rack for Bio-rad Evolis
- Made in Germany





RACKED FILTERED TIP

Filter tips protects the sample from aerosols and liquids. A proper barrier to all contaminants, the tips comes with reliable liquid blocking technology.



Available in: White Tips, 10µL Yellow Tips, 200µL Blue Tips, 1000 μL





- +632 8287.3683 +63917.878.2288 +63925.878.2288 www.trulaboratories.com