

ROTATOR

Digital Rotator DSR-2100D-N/DSR-2800D-N are widely used in VDRL, RPR and EIA test. Digital Rotator is equipped with digital display to indicate present speed and the remaining time.

FEATURES:

- Wide used in VDRL, RPR, and EIA test.
- Stable and quiet operation.
- Simple and Digital types for your options
- With variable speed, timed and continuous mode operation
- Quality Certification: ISO 9001, ISO 13485



SPECIFICATIONS

MODEL	DSR-2100P (ANALOG)	DSR-2100D-N (DIGITAL)
Speed	20~240 rpm	20~240 rpm
Time Mode	Continuous mode / timed mode 0~30 min. Mechanical timer)	Continuous mode / timed mode (99 min. 59 sec. Digital timer)
Display	NONE	Digital display for present speed value and remaining time
Platform Dimension	230 x 230 mm	230 x 230 mm
Orbit Diameter	20 mm	20 mm
Load Capacity	2.0 kgs. at at the speed from 20~150 rpm; 1.0 kgs. at the speed from 151~240 rpm	2.0 kgs. at at the speed from 20~150 rpm; 1.0 kgs. at the speed from 151~240 rpm
Accessory	1 pc. of spring	1 pc. of spring
Weight (Kg.) (N.W. / G.W.)	3.7 / 4.5	3.9 / 4.7
Overall Dimension (W x D x H)	260 x 295 x 140 mm	260 x 295 x 140 mm
Power Supply	AC 220V 60Hz 50/60Hz available	AC 220V 60Hz 50/60Hz available

PIPETTE SHAKER

The Pipette Shaker can easily achieve the uniform mixing of the chemicals with blood while they are in the diluting pipettes.

FEATURES:

- Has a pipette holder made of stainless steel spring clips.
- Provides continuous or timed operation (0-15 minutes)
- Can hold up to 6 pipettes
- Quality certificate: ISO 9001, ISO 13485



MODEL: PS-600T

SPECIFICATIONS:

MODEL	CAPACITY	SWITCH	FIXED SPEED	WEIGHT (KG) (N.W. / G.W.)	OVER ALL DIMENSIONS (W x D x H)
PS-600T	6 Pipettes	3 Positions (ON/ OFF / 0-15 min. timer)	2,500 rpm	1.8 / 2.2	150 x 135 x 160 mm

Power Supply: AC (220V 50/60Hz available)