

# Foresight<sup>®</sup> Thyroid EIA Test Kits TSH, Total T3, Free T3, Total T4 and Free T4



Foresight<sup>®</sup> Thyroid EIA Test Kits offer high quality and accurate Thyroid testing at an affordable price!

- *Accurate and Reliable*
- *Convenient*
- *Easy to Use*



# Foresight<sup>®</sup> Thyroid EIA Test Kits

## Accurate and Reliable

- High sensitivity with good reproducibility for all thyroid function tests
- Wide calibration range differentiates among various thyroid disorders
- Enhanced procedure available for TSH test as a helpful option for diagnosis of hyperthyroidism
- Cost-effective and comprehensive thyroid testing suitable for regular patient monitoring

## Convenient

- Standardized procedure for TSH, Total T3, Free T3, Total T4 and Free T4
- Consistent specimen volume and reagent volumes across all 5 tests
- Total incubation time as short as 75 minutes
- Room temperature incubation eliminates the need for incubators
- Ready to use reagents shorten initial preparation time

## Easy to Use

- Color-coded and numbered labeling for quick visual identification
- Individualized color for each reagent and across all EIA tests provides procedural control and efficient testing
- 96 microwell plate designed with breakaway 8-well strips in resealable foil pouch
- Compatible for use on automated sample processors or manual operation
- 12 month shelf life



Product Name	Catalog No.	Test Type	Incubation Time	Sensitivity	Tests per Kit*
Thyroid	TSH <sup>✓</sup>	I231-3011	Quantitative Standard Procedure: 75 min (60/15 min) Enhanced Procedure: 135 min (120/15 min)	Standard Procedure: 0.078 µIU/mL Enhanced Procedure: 0.03 µIU/mL	96
	Total T3 <sup>▲</sup>	I231-3041	Quantitative 75 min (60/15 min)	TBD	96
	Free T3 <sup>▲</sup>	I231-3051	Quantitative 75 min (60/15 min)	TBD	96
	Total T4 <sup>▲</sup>	I231-3021	Quantitative 75 min (60/15 min)	TBD	96
	Free T4 <sup>▲</sup>	I231-3031	Quantitative 75 min (60/15 min)	TBD	96

\* 480 and 48 tests per kit available soon    ▲ Coming soon

✓ **CE Marked for sale in the European Community** 

**Other EIA products available:**

HIV, HBsAg, HBsAb, HBeAg, HBeAb, HBcAb, HCV, HEV IgM<sup>✓</sup>, Syphilis<sup>✓</sup>, ToRCH, *H. pylori* Antigen<sup>✓</sup>, AFP<sup>✓</sup> and CEA<sup>✓</sup>.

Upcoming EIA products including HIV Ag/Ab, *H. pylori* IgG and IgM.

**We also offer other rapid diagnostic and medical products for:**

Clinical Chemistry including Urinalysis, Blood Glucose Monitoring and Others.

**Contact us for worldwide distribution and customer manufacturing (OEM) opportunities**



ACON Laboratories, Inc., 4108 Sorrento Valley Boulevard, San Diego, CA 92121, U.S.A.

Tel: 1-858-875-8000 • Fax: 1-858-535-2035 • E-mail: [info@aconlabs.com](mailto:info@aconlabs.com)

Please visit our website for details: [www.aconlabs.com](http://www.aconlabs.com)

Authorized Distributor